




**Hemo bioscience™**  
The science in Blood Banking

READY-TO-USE QUALITY CONTROL FOR GEL CARDS



QUALITY CONTROL FOR C.A.T.  
COLUMN AGGLUTINATION TECHNOLOGY

**Hemo bioscience MQC-CAT** is the **ONLY READY-TO-USE** product currently available that is intended for quality control testing of routine blood bank reagents used in manual column agglutination technology (CAT) blood typing and antibody screening systems.

The kit includes positive and negative red cells and plasma, providing a cost-effective method for quality control testing of routine CAT testing reagents and cards.

Evaluated with Grifols DG Gel® 8 and Ortho Clinical Diagnostics ID-MTS™ Gel systems.

PRODUCT NAME	DESCRIPTION	PRODUCT CODE	PRODUCT SIZE
MQC-CAT	Blood bank QC for CAT methods	H9350	RBC-POS 1 x 4 mL RBC-NEG 1 x 4 mL PLASMA-POS 1 x 10 mL PLASMA-NEG 1 x 10 mL

FEATURES & BENEFITS

- Unique ready to use product
- 63 day expiry dating
- Vials with color coded caps – easy to use
- Comprehensive worksheet included
- Red cell antigen profile provided – perform quality control of Rh antisera
- Compatible with Grifols DG Gel® 8 and Ortho Clinical Diagnostics ID-MTS™ Gel systems

## HEMO BIOSCIENCE COMMITMENT TO SERVICE

### CUSTOMER SERVICE

**Hemo bioscience** offers you flexible ordering options by fax, phone or email. To assist you in inventory management, we can set up standing orders or blanket orders. Most products are available for immediate overnight shipping. We develop, manufacture and inventory products in our facility, located in the heart of Research Triangle Park, North Carolina.

### TECHNICAL SUPPORT

**Hemo bioscience** technical support is staffed by professional, qualified technologists who have donor center, immunohematology reference laboratory and transfusion laboratory experience; we understand your challenges in the laboratory. We have two fully operational immunohematology laboratories on site with qualified staff who can resolve product questions and guide you in deciding which products fit your needs.

### QUALITY POLICY

**Hemo bioscience** strives to provide our customers with products and customer service that meet or exceed their expectations. We have established a comprehensive Quality Management System which provides a framework for measuring and improving our performance based on FDA CFR 21 part 820 standards, and we are ISO 13485 accredited. We have systems and procedures in place to support us in our goal of ensuring our customer requirements are met.



633 Davis Dr. Ste. 100  
Research Triangle Park  
Morrisville, NC 27560

Tel: 1 866 332 2835 info@hemobioscience.com  
Tel: 1 919 313 2888 www.hemobioscience.com  
Fax: 1 919 313 2894

