



**Hemo bioscience**<sup>™</sup>  
*The science in Blood Banking*

## ELUclear-ELUTION KIT

**ELUclear**, from **Hemo bioscience**, is a glycine acid red cell elution kit that is used to elute IgG antibodies from red cell membranes so that they can be identified. The ELUclear kit contains an eluting solution with a pH indicator, a buffering solution, and a concentrated wash solution which minimizes dissociation during the wash phase.

The eluting solution and buffering solution are supplied in color-coded dropper vials for ease of use and to help prevent cross contamination. For added convenience, a wash bottle for working wash solution preparation is included with every kit.

Eluates prepared from the ELUclear kit can be tested using both traditional tube and column agglutination technology (CAT) testing methods.

PRODUCT NAME	DESCRIPTION	PRODUCT CODE	PRODUCT SIZE
ELUclear	Glycine acid elution kit for red cells	H360	Sufficient for 10 x 1 mL elutions

### FEATURES & BENEFITS

- Clearly defined endpoint color change
- Validated for both tube and gel techniques
- Acid and buffer provided in dropper vials to reduce risk of contamination
- Color-coded caps help to differentiate step 1 (acid) and step 2 (buffer)
- Wash bottle with fill line allows for easy preparation of working wash solution

## HEMO BIOSCIENCE COMMITMENT TO SERVICE

### CUSTOMER SERVICE

**Hemo bioscience** offers you flexible ordering options by fax, phone or email. To assist you in inventory management, we can set up standing orders or blanket orders. Most products are available for immediate overnight shipping. We develop, manufacture and inventory products in our facility, located in the heart of Research Triangle Park, North Carolina.

### TECHNICAL SUPPORT

**Hemo bioscience** technical support is staffed by professional, qualified technologists who have donor center, immunohematology reference laboratory and transfusion laboratory experience; we understand your challenges in the laboratory. We have two fully operational immunohematology laboratories on site with qualified staff who can resolve product questions and guide you in deciding which products fit your needs.

### QUALITY POLICY

**Hemo bioscience** strives to provide our customers with products and customer service that meet or exceed their expectations. We have established a comprehensive Quality Management System which provides a framework for measuring and improving our performance based on FDA CFR 21 part 820 standards, and we are ISO 13485 accredited. We have systems and procedures in place to support us in our goal of ensuring our customer requirements are met.



633 Davis Dr. Ste. 100  
Research Triangle Park  
Morrisville, NC 27560

Tel: 1 866 332 2835 info@hemobioscience.com  
Tel: 1 919 313 2888 www.hemobioscience.com  
Fax: 1 919 313 2894

