


HEMO-QC

DAILY QUALITY CONTROL FOR YOUR BLOOD BANK

Hemo-QC is a whole blood quality control kit and is intended for the daily quality control of routine blood bank reagents used in manual, semi-automated and automated blood typing and antibody screening test systems.

The kit consists of 3 tubes that provide a cost-effective method of daily QC and confirm the correct performance of ABO, Rh(D) and antibody detection reagents.



PRODUCT NAME	DESCRIPTION	PRODUCT CODE	PRODUCT SIZE	MONTHLY USE
Hemo-QC	Quality control for blood banking reagents	H9200	1 kit = 3 x 7 mL tubes	3-4 kits/month

FEATURES & BENEFITS

- 7 mL tubes are compatible with all instrument systems currently available in the United States
- 3 tubes perform daily QC
- 42-day closed vial stability and 10-day open vial stability
- Cost effective – fewer reagent sets used monthly
- Monthly shipments on standing order

HEMO BIOSCIENCE COMMITMENT TO SERVICE

CUSTOMER SERVICE

Hemo bioscience offers you flexible ordering options by fax, phone or email. To assist you in inventory management, we can set up standing orders or blanket orders. Most products are available for immediate overnight shipping. We develop, manufacture and inventory products in our facility, located in the heart of Research Triangle Park, North Carolina.

TECHNICAL SUPPORT

Hemo bioscience technical support is staffed by professional, qualified technologists who have donor center, immunohematology reference laboratory and transfusion laboratory experience; we understand your challenges in the laboratory. We have two fully operational immunohematology laboratories on site with qualified staff who can resolve product questions and guide you in deciding which products fit your needs.

QUALITY POLICY

Hemo bioscience strives to provide our customers with products and customer service that meet or exceed their expectations. We have established a comprehensive Quality Management System which provides a framework for measuring and improving our performance based on FDA CFR 21 part 820 standards, and we are ISO 13485 accredited. We have systems and procedures in place to support us in our goal of ensuring our customer requirements are met.



633 Davis Dr. Ste. 100
Research Triangle Park
Morrisville, NC 27560

Tel: 1 866 332 2835 info@hemobioscience.com
Tel: 1 919 313 2888 www.hemobioscience.com
Fax: 1 919 313 2894

