



Hemo bioscience™
The science in Blood Banking

BLOOD BANK MQC (MANUAL QUALITY CONTROL)



REAGENT QUALITY CONTROL

Hemo bioscience MQC is intended for quality control testing of routine blood bank reagents used in manual blood typing and antibody screening systems.

The kit includes positive and negative red cells and plasma, providing a cost-effective method for quality control testing of routine reagent red blood cells and antisera.

PRODUCT NAME	DESCRIPTION	PRODUCT CODE	PRODUCT SIZE
MQC	Blood bank QC for manual methods	H9300	RBC-POS 1 x 10 mL RBC-NEG 1 x 10 mL PLASMA-POS 1 x 10 mL PLASMA-NEG 1 x 10 mL

FEATURES & BENEFITS

- Cost effective – 63 day expiry
- Dropper vials with color coded caps – prevent cross contamination
- Worksheet included for recording results
- Can be used to perform quality control testing of rare antisera (C,c,E,e,K,Fy^a,Fy^b,Jk^a,Jk^b) – antigen profile provided for red cells

HEMO BIOSCIENCE COMMITMENT TO SERVICE

CUSTOMER SERVICE

Hemo bioscience offers you flexible ordering options by fax, phone or email. To assist you in inventory management, we can set up standing orders or blanket orders. Most products are available for immediate overnight shipping. We develop, manufacture and inventory products in our facility, located in the heart of Research Triangle Park, North Carolina.

TECHNICAL SUPPORT

Hemo bioscience technical support is staffed by professional, qualified technologists who have donor center, immunohematology reference laboratory and transfusion laboratory experience; we understand your challenges in the laboratory. We have two fully operational immunohematology laboratories on site with qualified staff who can resolve product questions and guide you in deciding which products fit your needs.

QUALITY POLICY

Hemo bioscience strives to provide our customers with products and customer service that meet or exceed their expectations. We have established a comprehensive Quality Management System which provides a framework for measuring and improving our performance based on FDA CFR 21 part 820 standards, and we are ISO 13485 accredited. We have systems and procedures in place to support us in our goal of ensuring our customer requirements are met.



633 Davis Dr. Ste. 100
Research Triangle Park
Morrisville, NC 27560

Tel: 1 866 332 2835 info@hemobioscience.com
Tel: 1 919 313 2888 www.hemobioscience.com
Fax: 1 919 313 2894

