



Hemo bioscience™
The science in Blood Banking

MEDTEK – MEDICAL TECHNOLOGY EDUCATION KIT

RESOURCES FOR THE STUDENT LABORATORY

Simulated Patient Plasma (SPP) from **Hemo bioscience** was designed with the student or training laboratory in mind. MedTEK SPP is a cost-effective, ready-to-use product that accurately simulates a patient sample containing clinically significant IgG alloantibodies. Standardized for use with all major testing techniques, these human sourced materials are an effective and convenient teaching resource.

PRODUCT NAME	SPECIFICITIES	VOLUME
MedTEK SPP	anti-D, C, c, E, K, S, Fya, Fyb	75 ml

*Materials found non-reactive for HBsAg, anti-HIV 1/2, and anti-HCV.
 These products are not for in vitro diagnostic use and are not FDA cleared.*

FEATURES & BENEFITS

- IgG antibodies show 2-3+ reactivity
- Suitable for tube and gel AGT techniques
- SPP specificities can be blended
- Customizable ABO type; anti-A and anti-B concentrate provided with each order

EDUCATIONAL SAVINGS!

We now offer a discount of **15% off list price** on all our products to institutions of higher learning that offer MT, MLT and CLS certifications. Call us for details on how to qualify.

(866) 332-2835

HEMO BIOSCIENCE COMMITMENT TO SERVICE

CUSTOMER SERVICE

Hemo bioscience offers you flexible ordering options by fax, phone or email. To assist you in inventory management, we can set up standing orders or blanket orders. Most products are available for immediate overnight shipping. We develop, manufacture and inventory products in our facility, located in the heart of Research Triangle Park, North Carolina.

TECHNICAL SUPPORT

Hemo bioscience technical support is staffed by professional, qualified technologists who have donor center, immunohematology reference laboratory and transfusion laboratory experience; we understand your challenges in the laboratory. We have two fully operational immunohematology laboratories on site with qualified staff who can resolve product questions and guide you in deciding which products fit your needs.

QUALITY POLICY

Hemo bioscience strives to provide our customers with products and customer service that meet or exceed their expectations. We have established a comprehensive Quality Management System which provides a framework for measuring and improving our performance based on FDA CFR 21 part 820 standards, and we are ISO 13485 accredited. We have systems and procedures in place to support us in our goal of ensuring our customer requirements are met.



633 Davis Dr. Ste. 100
Research Triangle Park
Morrisville, NC 27560

Tel: 1 866 332 2835 info@hemobioscience.com
Tel: 1 919 313 2888 www.hemobioscience.com
Fax: 1 919 313 2894

